



art.-no.: 75-091-100

type: **Cathodic Disbonding Tester**
multifunctional tester for 4 measuring stations

acc. to current standards: ASTM G 8 – ASTM G 42 – ASTM G 80 – ASTM G 95
ISO EN 15 711
DIN 30 671

description: multi-channel potentiostat/galvanostat for control up to 4 measuring stations as well as voltage and power controlled
LCD display for current max. 400 mA, 3 range, switchable display 400 mA/40 mA/4 mA
LCD display for potential/voltage
voltage range adjustable ± 15 volt
data acquisition of the measuring values via SUB-D-system plug possible
set of 4 connection cables with plug

technical data: weight: approx. 14 kg
dimensions (WxHxD): 400 x 240 x 280 mm
power supply: 230 V, 50/60 Hz

