



AUTOMATIC CLOUD AND POUR POINTS

ASTM D97 - ASTM D2500 - NF EN 23015 - ISO 3015 - NFT 60105

These testers allow to test automatically Cloud and Pour Points and both connect up by means of a mixed apparatus.

Fully in accordance with above standards.

Sample tipping during the Pour Point measure and maintenance of the test tube in horizontal position for 5 seconds to validate the end of test.

CHARACTERISTICS

- Automatic piloting and management by a process control unit
- Programmable pre-heating of sample
- Easy use and maintenance
- Standardized method and fast research

SPECIFICATION

The Automatic Cloud and Pour Points may easily be integrated into a Quality Control Procedure. The analytic concept is the same as ONLINE analysers, and is strictly conform to the standards.