



OIL SEPARATION FROM LUBRICATING GREASE DURING STORAGE

ASTM D1742

DEFINITION

This test method covers the determination of the tendency of lubricating grease to separate oil during storage

CHARACTERISTICS

- Temperature maintenance by electronic regulation
- Stiff manifold
- Air pressure adjustment
- Pressure gauge

SPECIFICATION

6 - unit assembled Tester including constant temperature oven, 6-test assemblies A type with cup, cover, 200 mesh sieve, glass beaker, manifold with manometer and reducing valve, stopcocks and ASTM 5C thermometer