



VISCOMETERS

SAYBOLT ASTM D 88 - ASTM D 244 - ASTM E 102

ENGLER ASTM D 1665 - DIN 51560

REDWOOD IP 70

DEFINITION

Those testers allow the determination of Saybolt, Engler or Redwood viscosities. They do not determine absolute viscosity of petroleum products, tar or fluid products, but an empirical flow test, only by conforming to requirements of the test methods allowing results which can be reproducible.

CHARACTERISTICS

- Polyvalent apparatus
- Electric heating
- Automatic control of temperature
- Stirring motor for water circulation

SPECIFICATION

Those testers are equipped with electric heater, digital display of temperature, stirring motor for bath temperature standardization, possibility to cool down the bath by means of an external cooling system.