



## **MANUAL TAG CLOSED CUP TESTER**

### **ASTM D56, IP 304**

This test method covers determination of flash point by TAG closed tester of liquids with a viscosity below 9,5 cst at 25°C or 5,5 cst at 40°C, and a flash point below 93°C, except cut-back asphalts.

#### **CHARACTERISTICS**

- Electric heating
- Temperature controlled bath
- Adjustment of heating power
- Gas ignition device

#### **SPECIFICATION**

MANUAL TAG CLOSED CUP TESTER comprises a complete test assembly, a cup, a cover, a pilot light and a test flame with gas ignition, an electric heating with regulator, a bath and ASTM 9C and 57°C thermometers.