


# RAY-RAN

# POLYTEST

## Advanced Microprocessor Controlled Density Column Filling System



- Extremely Quick and accurate compared with conventional systems.
- Any two miscible liquids can be used to give lighter or heavier densities.
- No Pre-adjustment of density is necessary for the top and bottom of the column.
- Column volume and upper / lower densities are fully programmable.
- Variable buffer volumes can be programmed at the top and bottom of the column.
- Details of previous column settings are automatically stored for re-use.
- 4 speed column sweep mechanism facility.
- Software automatically calculates density of objects once column has been calibrated.



**RAY-RAN**  
TEST EQUIPMENT LTD

WORLD WIDE TESTING FOR QUALITY